

Field-proven CAD Technology

Optimal Clinical Flexibility

ImageChecker[®] 2D
Computer-Aided Detection

We understand the constantly changing digital environment and the need for a flexible, compatible CAD solution. ImageChecker[®] delivers proven performance and has been the leading CAD technology for more than a decade.

Using sophisticated software algorithms, ImageChecker searches digital mammograms for potential microcalcifications and masses—characteristics commonly associated with breast cancer. Marks are placed on the regions-of-interest to focus the radiologist's attention. The highly-evolved algorithms, refined to deliver extremely sensitive results without excessive false-positive marks, provide streamlined case review.

Advanced Technology

In addition to detecting regions-of-interest, the ImageChecker algorithms also incorporate anatomic correlation technology. This advanced technology analyzes corresponding findings in the contralateral breast and different views of the same breast.

Flexible Solutions

ImageChecker provides optimal flexibility by allowing you to select sensitivity settings, or operating points, that are most suitable for your needs. ImageChecker ranks findings in order of the prominence of features in a given region. A region will be marked only if the ranking falls above a chosen operating point. Clinical practices can select from three different operating points each for calcification and mass detection, for a total of nine combinations.

ImageChecker supports images from the following digital mammography systems:

| Manufacturer | USA | Other |
|--------------------|-----|-------|
| Hologic | Y | Y |
| GE Medical Systems | Y | Y |
| Siemens | Y | Y |

ImageChecker operates on the Cenova™ server, providing complete compatibility with DICOM conformant networks.


Field-proven CAD Technology

Optimal Clinical Flexibility




ImageChecker® 9.3 Algorithm Performance¹

The following graphs show plots of algorithm sensitivity vs false-mark rate based on cases with the four standard screening views, with data points for each of the three operating points.



Calcification Performance



Mass Performance

An operating point of “0” corresponds to the lowest sensitivity and the fewest false-marks; an operating point of “2” corresponds to the highest sensitivity and the most false-marks. Practices may select separate operating points for calcifications and masses.

Operating Points

| | Calcifications | | | Masses | | |
|--------------------------|----------------|------|----------------|--------|----------------|-----|
| | 0 | 1 | 2 | 0 | 1 | 2 |
| Sensitivity ² | 95% | 96% | 97% | 83% | 88% | 90% |
| False Marks ³ | 0.36 | 0.48 | 0.64 | 0.68 | 1.0 | 1.4 |
| | | | <i>Default</i> | | <i>Default</i> | |

Operating Point Combinations

| Mass/Calc | 0/0 | 1/1 | 2/2 |
|-------------------------------|-----|-----|-----|
| Case specificity ⁴ | 48% | 35% | 24% |

Operating Environment

Operating Platform: Cenova™
 Output Supported: PACS or DICOM workstations

Performance on Cenova

Maximum Number of Digital Ports: 4 (optional)
 Case Throughput (four-image): 30-60 cases/hr (nominal)⁶

Extended Features

RightOn™ CAD Marks
 Citra® Core⁵

Assorted “shaped” markers indicate the types of features that were detected

- EmphaSize™ marks are scaled according to the prominence of features
- PeerView® Digital provides anatomic outlines of tissue
- LesionMetrics™ provide region specific information, such as lesion *Size, Distance to chest wall and Distance to nipple*

¹ Performance derived from analog data. ImageChecker CAD algorithms have been assessed using both digitized film mammograms and digital mammograms. Algorithm performance between the two has been shown to be comparable.

² Data based on 1355 biopsy proven breast cancers (767 mass and 588 calcification cases)

³ Data based on 445 four-view normal cases

⁴ Four film normal cases with no markers

⁵ Available only on advanced SecurView workstations.

⁶ Performance dependent on recommended hardware, network bandwidth and input rate of images

United States
 35 Crosby Drive
 Bedford, MA 01730-1401 USA
 Tel: +1.781.999.7300
 Sales: +1.781.999.7453
 Fax: +1.781.280.0668
 www.hologic.com

To learn more about ImageChecker please contact your Hologic account manager, call **781.999.7300** or email **womenshealth@hologic.com**



DS-00025 US (10/10) ©2010 Hologic, Inc. All Rights Reserved. Printed in the United States. Specifications are subject to change without prior notice. Hologic, ImageChecker, RightOn, EmphaSize, PeerView, LesionMetrics, Citra, and associated logos are trademarks and/or registered trademarks of Hologic, Inc. and/or its subsidiaries in the United States and other countries.