



**Affirm™**  
Breast Biopsy Guidance System



Affirm™  
Breast Biopsy Guidance System

## Advanced Solutions in Upright Breast Biopsy

The Affirm™ breast biopsy guidance system represents the next generation in upright breast biopsy. This innovative solution is designed for our Selenia® Dimensions® digital mammography system to streamline efficiency and optimize workflow, and gives facilities the opportunity to offer leading-edge technology in minimally invasive breast biopsy. With a quick and easy transition from mammography to upright biopsy, the Affirm system offers practices a cost-effective, space-saving and ergonomic system to expand its range of service offerings.



Selenia Dimensions digital mammography system with Affirm breast biopsy guidance system

### A Comprehensive Solution

Hologic® is the first to deliver a comprehensive breast health solution featuring:

- Direct capture digital mammography.
- The world's first commercially available breast tomosynthesis system.
- Upright and prone breast biopsy guidance systems.
- Vacuum-assisted breast biopsy devices utilizing the widest variety of needle sizes.
- Hologic biopsy site markers.

**Hologic breast health solutions allow physicians to meet the needs of the broadest spectrum of patients.**



ATEC® and Eviva® breast biopsy devices



SecurMark® biopsy site markers



## Innovative Biopsy Solution

The Affirm™ breast biopsy guidance system brings new levels of performance to upright interventional procedures. Its novel, 10° angled biopsy approach provides an unobstructed view of lesions. With other innovations including automated image acquisition, sophisticated yet easy-to-use targeting software and advanced ergonomic design, the Affirm system is poised to expand the horizons of breast biopsy.

## Streamlined Workflow

- **Automated Image Acquisition**

Only the Selenia® Dimensions® digital mammography system with the Affirm system offers one easy step to automate stereo image acquisition: When the exposure switch is pressed, the x-ray tube automatically moves to the correct position to expedite image acquisition and reduce procedure time.

- **Fully integrated user interface**

All activities are performed on the Selenia Dimensions acquisition workstation, providing easy to use, intuitive controls and simplified workflow.

- **Accurate and efficient targeting**

Cartesian targeting software removes guesswork and provides visual feedback of needle placement.





5



### **Flexibility and Compatibility**

The Affirm™ system's versatile design enables use with a wide range of biopsy devices. Also, its pre-programmed settings for Hologic® Eviva® and ATEC® biopsy devices offer quick and efficient set-up for reduced procedure times.

The Affirm system simplifies patient positioning with an adjustable breast biopsy imaging field of view (up to 18 cm x 24 cm). Our advanced detector provides superb image quality.

The Affirm system's paddles are designed with biopsy in mind. Crystal clear, smooth-edged paddles offer greater patient comfort, enhanced breast tissue grip, and a full field of view. The paddle kit includes standard biopsy paddles and a specially designed paddle for managing difficult axillary positioning.

These exclusive features ensure that the Selenia® Dimensions® digital mammography system with the Affirm system will meet all your patient needs.

### **Ergonomic Design**

Weighing approximately 7 kg (15 lb), the Affirm system's compact design and ergonomic, easy-grip handles let you install the system in just seconds. Easily transportable, the Affirm system is especially convenient for centers using multiple Selenia Dimensions digital mammography systems, as it can be licensed for use with more than one gantry.

Its integrated design leverages the 70 cm SID to provide a comfortable working space and larger clearance for installation of biopsy devices.

The biopsy control module can be positioned on either side of the Affirm system, ensuring optimal access. Touch screen controls and simplified user interface offer added efficiency during procedures.

### **A Revolution in Breast Biopsy**

The Affirm system is a next-generation biopsy guidance system that allows you to perform interventional procedures under the same modality as diagnostic procedures.

The Selenia Dimensions digital mammography system with the Affirm system, the industry's most advanced breast biopsy guidance system, lets you offer your patients a comprehensive breast care solution.

## Optimal Access and Visualization.

The Affirm™ breast biopsy guidance system was developed to address the challenges of upright breast biopsy. This novel solution was designed to meet the biopsy needs of today and paves the way for future advances in interventional procedures. The Affirm system delivers:

- A fully integrated solution optimized for the Selenia® Dimensions® digital mammography system that minimizes procedure steps and simplifies workflow.
- Intuitive user interface driven by the Selenia Dimensions acquisition workstation for enhanced ease of use.
- Advanced ergonomic and lightweight design for quick and easy transition from diagnostic to interventional procedures under the same imaging modality.
- 70 cm SID, the longest in the industry, for easier positioning of patients and installation of biopsy devices.
- A unique, 10° angled biopsy approach for an unobstructed view of lesions while maintaining simple Cartesian targeting methods in all biopsy modes.
- Compatibility with a wide array of biopsy devices for greater convenience.
- A complete solution for screening, diagnostic, and interventional procedures.



**Hologic is defining the standard of care in women's health. Our technologies help doctors see better, know sooner, reach further and touch more lives. At Hologic, we turn passion into action, and action into change.**

BREAST IMAGING SOLUTIONS | INTERVENTIONAL BREAST SOLUTIONS | BONE HEALTH  
PRENATAL HEALTH | GYNECOLOGICAL HEALTH | MOLECULAR DIAGNOSTICS

**HOLOGIC®**

The Women's Health Company

**United States / Latin America**

35 Crosby Drive  
Bedford, MA 01730-1401 USA  
Tel: +1.781.999.7300  
Sales: +1.781.999.7453  
Fax: +1.781.280.0668  
www.hologic.com

**Europe**

Everest (Cross Point)  
Leuvensesteenweg 250A  
1800 Vilvoorde, Belgium  
Tel: +32.2.711.4680  
Fax: +32.2.725.2087

**Asia Pacific**

7th Floor, Biotech Centre 2  
No. 11 Science Park West Avenue  
Hong Kong Science Park  
Shatin, New Territories  
Hong Kong  
Tel: +852.3748.7700  
Fax: +852.3526.0723

**Australia**

Hologic (Australia) Pty Ltd  
Suite 402, Level 4  
2 Lyon Park Road  
Macquarie Park NSW 2113  
Australia  
Tel: +61.2.9888.8000  
Fax: +61.2.9870.7555